



REF CS-12123-F

12 Fr. 3 Lumen 16 cm catheter length

.035 inch dia. spring-wire guide

## Large-Bore Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions

### Contents:

- 1: Multi-Lumen Indwelling Catheter: 12 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 3: Dust Cap: Non-Vented
- 1: Drape: 24" x 36" (61 cm x 91 cm)
- 1: Scalpel: #11

- fr** Jeu de cathétérisme veineux central de grand calibre à lumières multiples avec cathéter Blue FlexTip pour les perfusions de grand volume
- de** Besteck für viellumigen, zentralen, großkalibrigen Venenkatheter mit Blue-FlexTip-Katheter für Infusionen großer Flüssigkeitsmengen
- it** Set per cateterismo della vena centrale multilume di diametro ampio con catetere Blue FlexTip per volumi di infusione elevati
- pl** Wielokanałowy zestaw o dużej średnicy do cewnikowania żył centralnych z cewnikiem Blue FlexTip do przetaczania dużych objętości płynów infuzyjnych
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Grande Calibre, com Cateter Blue FlexTip, para Infusões de Grandes Volumes
- ru** Комплект для катетеризации главных вен с многопросветным катетером большого внутреннего диаметра с Blue FlexTip для вливаний большого объема
- sl** Komplet za centralno venosko katetizacijo z velikimi luknjami in več svetlinami s katetrom Blue FlexTip za infundiranje velike prostornine
- es** Conjunto de cateterización venosa central de luz múltiple y gran diámetro con catéter Blue FlexTip para infusiones de gran volumen
- sv** Set för storhålig flerkanalig central venkatetrering med Blue FlexTip-kateter för infusioner med hög volym
- tr** Büyük Açıklıklı Multi-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile, Yüksek Hacimli infüzyonlar için

All components are CE 0086 unless otherwise noted.


**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

### California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.  
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



\* Priming volumes are approximate and are done without accessories.  
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
www.teleflex.com

LBL028089 (02/20)



REF CS-12123-F

12 Fr. 3 Lumen 16 cm catheter length .035 inch dia. spring-wire guide

Large-Bore Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions