



REF CS-17702-E

7 Fr. 2 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

Rx only. **Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.


California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm. www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- fr** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip
- de** Besteck für zweilumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter
- it** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip
- pl** Dwukanałowy zestaw do cewnikowania żył centralnych z cewnikami Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- ru** Комплект для катетеризации главных вен с двухпросветным катетером Blue FlexTip
- sl** Komplet za centralno vensko katetrizacijo z dvema svetlinama s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip
- sv** Set för tvåkanalig central venkatetrering med Blue FlexTip-kateter
- tr** İki Lümenli Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (14 Ga.)	0.57	5000
Proximal (18 Ga.)	0.40	1500

* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



LOT



CE REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

LBL028165 (02/20)



REF CS-17702-E

7 Fr. 2 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter