

ARROW®



REF EU-14553-N (IPNXXXXXX)

5.5 Fr. 3 Lumen 30 cm catheter length .018 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter for Femoral Vein Insertion


Contents:

- 1: Multi-Lumen Indwelling Catheter: 5.5 Fr. x 12" (30 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
- 1: Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Safety Scalpel: #11

All components are CE 0086 unless otherwise noted.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (20 Ga.)	0.38	800
Medial (22 Ga.)	0.35	300
Proximal (22 Ga.)	0.39	300

* Priming volumes are approximate and are done without the injection cap.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

- be Набор з шматпросветным цэнтральным веназным катэтарам з наканечнікам Blue FlexTip для катэтарызацыі сцёгнавай вены
- hr Komplet za centralnu vensku kateterizaciju s više lumena i kateterom Blue FlexTip za uvođenje u femoralnu venu
- cs Vicelumenový centrální žilní katetrizační set s katetrem Blue FlexTip, pro zavedení do femorální žíly
- da Multilumen centralvenekateteriseringssæt med Blue FlexTip kateter til indføring i v. femoralis
- nl Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip katheter voor r aanbrengen in de vena femoralis
- et Reieveeni sisestamiseks mõeldud mitme valendikuga tsentraalveeni kateteriseerimise komplekt sinise FlexTip kateetriga
- fi Suuraukkoinen monilumeninen keskuskaskimokatetri sarja, jossa on FlexTip katetri reisilaskimoon asettamiseen
- fr Set de cathétérisme veineux central multi-lumières avec cathéter Blue FlexTip pour une insertion par la veine fémorale
- de Zentralvenenkatheter-Set mit Multilumen mit Blue FlexTip Katheter zur Einführung in die V. femoralis
- el Σετ Καθετηριασμού Κεντρικής Φλέβας Πολλαπλού Αυλού με Καθετήρα Blue FlexTip για Εισαγωγή στη Μηριαία Φλέβα
- hu Többlumenű, centrális vénás katéterezési szett, Blue FlexTip káteeterrel a vena femoralis felőli bevezetéshez
- is Fjölrása miðlægt holæðaleggssett með Blue FlexTip æðalegg til isetningar í lærleggsbláæð
- it Set per cateterismo della vena centrale multilume con catetere Blue FlexTip per inserimento nella vena femorale
- lt Daugiakanalis centrinės venos kateterizavimo rinkinys su „Blue FlexTip“ kateteriu, skirtas įvesti į šlaunies veną
- no Multilumen Sentralt venekateteriseringssett med Blue FlexTip kateter for innføring i lårbensvene
- pl Zestaw do cewnikowania żył centralnych z wielokanalowym cewnikiem z końcówką Blue FlexTip, do wprowadzenia przez żyłę udową
- pt Conjunto de Cateterização Venoso Central com Lúmen Múltiplo e Cateter Blue FlexTip para Inserção na Veia Femoral
- ro Set de cateterizare venoasă centrală multi-lumen cu cateter cu Blue FlexTip pentru introducere în vena femurală
- ru Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip для катетеризации бедренной вены
- sr Komplet sa više lumena za centralnu vensku kateterizaciju sa plavim FlexTip kateterom za uvođenje u femoralnu venu
- sk Viačlúmenová centrálna venózna katetrizačná súprava s katétrom Blue FlexTip na zavedenie do stehennej žily
- sl Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip za vstavljanje v femoralno veno
- es Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip para introducción en la vena femoral
- sv Centralvenös katetersats med flera lumen och blå FlexTip kateter införing i vena femoralis
- tr Femoral Ven Girişi için Blue FlexTip Kateter ile Çok Lümenli Santral Venöz Kateterizasyon Seti
- uk Набір з багатопросвітним центральним венозним катетером з кінчиком Blue FlexTip для катетеризації стегнової вени



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

CE REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

LBL028168 (02/20)



REF EU-14553-N 5.5 Fr. 3 Lumen 30 cm catheter length .018 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter for Femoral Vein Insertion

ARROW®