



**REF** EU-16402-N (IPNXXXXXX)

**4** Fr. **2** Lumen **30** cm catheter length **.018** inch dia. spring-wire guide

## Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter for Femoral Vein Insertion

**Contents:**


- 1: Two-Lumen Indwelling Catheter: 4 Fr. x 12" (30 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
- 1: Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Safety Scalpel: #11

All components are CE 0086 unless otherwise noted.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- be** Набор з дзьвухпросветным цэнтральным вьянозным катэтарам с наканечнікам Blue FlexTip для катэтарызацыі сцягнавой вены
- hr** Komplet za centralnu vensku kateterizaciju s dva lumena i kateterom Blue FlexTip za uvođenje u femoralnu venu
- cs** Dvoulumenový centrální žilní katetrizační set s katetrem Blue FlexTip, pro zavedení do femorální žily
- da** To-lumen centralvenekateteriseringssæt med og Blue FlexTip kateter indføring i v. femoralis
- nl** Twee-lumen centraal veneuze katheterisatieset met Blue FlexTip katheter voor aanbrengen in de vena femoralis
- et** Reieveeni sisestamiseks mõeldud kahe valendikuga tsentraalveeni kateteriseerimise komplekt sinise FlexTip katetriga
- fi** Kaksiaukkoinen keskuskaskimokatetri sarja, jossa on sininen FlexTip katetri reisilaskimoon asettamiseen
- fr** Set de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip pour une insertion par la veine fémorale
- de** Zentralvenenkatheter-Set mit Doppel-Lumen mit Blue FlexTip Katheter zur Einführung in die V. femoralis
- el** Σετ Καθετηριασμού Κεντρικής Φλέβας Διπλού Αυλού με Καθετήρα Blue FlexTip για Εισαγωγή στη Μηριαία Φλέβα
- hu** Kétlumenes központi vénás katéterkészlet Blue FlexTip katéterrel a vena femoralis felőli bevezetéshez
- is** Tveggja rása miðlæggt holæðaleggsett með Blue FlexTip æðalegg til isetningar í lærleggssbláæð
- it** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip per inserimento nella vena femorale
- lt** Dvieju kanalų centrinės venos kateterizavimo rinkinys su „Blue FlexTip“ kateteriu, skirtas įvesti į šlaunies veną
- no** To-Lumen innsetningssett for sentralt venekateter med Blue FlexTip kateter for innsetting i lårbensvene
- pl** Zestaw do cewnikowania żył centralnych z dwukanałowymi cewnikami z końcówką Blue FlexTip, do wprowadzenia przez żyłę udową
- pt** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip para Inserção na Veia Femoral
- ro** Set de cateterizare venoasă centrală cu două lumene cu cateter cu Blue FlexTip pentru introducerea în vena femorală
- ru** Набор с двухпросветным центральным венозным катетером с кончиком Blue FlexTip для катетеризации бедренной вены
- sr** Komplet sa dva lumena za centralnu vensku kateterizaciju sa plavim FlexTip kateterom za uvođenje u femoralnu venu
- sk** Dvojlúmenová centrálna venózna katetrizačná súprava s katétrom Blue FlexTip na zavedenie do stehennej žily
- sl** Dvolumenski centralni venski pribor za katetrizacijo s katetrom Blue FlexTip za vstavljanje femoralno veno
- es** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip para introducción en la vena femoral
- sv** Centralvenös katetersats med två lumen och blå FlexTip kateter för införing i vena femoralis
- tr** Femoral Ven Girişi için Blue FlexTip Kateter ile İki Lümenli Santral Venöz Kateterizasyon Seti
- uk** Набір з двохпросвітним центральним венозним катетером з кінчиком Blue FlexTip для катетеризації стегнової вени



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (22 Ga.)	0.30	190
Proximal (22 Ga.)	0.38	280

\* Priming volumes are approximate and are done without the injection cap.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



**LOT** LotNumber  
YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber

**CE** 0086  
Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
www.teleflex.com  
LBL048345 (02/20)



**REF** EU-16402-N

**4** Fr. **2** Lumen **30** cm catheter length **.018** inch dia. spring-wire guide

Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter for Femoral Vein Insertion

