



REF AK-17702-J

ARROW®

7 Fr. 2 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Two-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter

Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 22 Ga. x 2-3/4" (7 cm) TW
- 1: Syringe: 3 mL Luer-Lock
- 1: Syringe: 5 mL Luer-Slip
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 2: UserGard® Needle-Free Injection Hub
- 2: Dust Cap: Non-Vented
- 2: Injection Site Cap
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway® Disposal Cup
- 1: Drape: 21" x 36" (53 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive
- 1: Scalpel: #11
- 2: Stopcock: 3-way
- 2: Prep Sponge Swab
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle
- 1: Tray: Prep

All components are CE 0086 unless otherwise noted.


Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- fr** Kit de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip
- de** Ausstattung für zweilumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter
- it** Kit per cateterismo della vena centrale a due lumi con catetere Blue FlexTip
- pt** Kit de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- es** Juego de cateterización venosa central de dos luces con catéter Blue FlexTip
- sv** Uppsättning för tvåkanalig central venkatetrering med Blue FlexTip-kateter



Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA



Teleflex®
 www.teleflex.com
 LBL048356 (03/20)



REF AK-17702-J

7 Fr. 2 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Two-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter

ARROW®