



REF CS-22855 (IPNXXXXXX)

8.5 Fr. 5 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Arrowgard Blue® Five-Lumen CVC

Contents:

- 1: Arrowgard Blue® Five-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 16 cm
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 5: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Scalpel: #11

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

Not made with natural rubber latex.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F). Fluid path components are non-pyrogenic.

Contraindications: The Arrowgard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.

All components are CE 2797 unless otherwise noted.



CAS # 7440-48-4



| Lumen | Priming Volume* (mL) | Gravity Flow Rate† (mL/hr) |
|-------------------|----------------------|----------------------------|
| Distal (16 Ga.) | 0.58 | 3030 |
| Medial 1 (14 Ga.) | 0.71 | 5470 |
| Medial 2 (18 Ga.) | 0.36 | 1046 |
| Medial 3 (18 Ga.) | 0.38 | 1187 |
| Proximal (18 Ga.) | 0.40 | 1212 |

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

EU Authorized Representative and Importer:

Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of xxx
Packaged in xxx

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Teleflex®
www.teleflex.com

LBL051195 R02 (2025-01)



8.5 Fr. 5 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

REF CS-22855

Arrowgard Blue® Five-Lumen CVC

