

Product of xxx
Packaged in xxx



REF BR-22703 (IPNXXXXXX)

7 Fr. 3 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Arrowgard Blue® Three-Lumen CVC

Contents:

- 1: Arrowgard Blue® Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 16 cm
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Scalpel: #11
- 3: Stopcock: 4-way High-Flow

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

Not made with natural rubber latex.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F). Fluid path components are non-pyrogenic.

Contraindications: The Arrowgard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.

All components are CE 2797 unless otherwise noted.



CAS # 7440-48-4



| Lumen | Priming Volume* (mL) | Gravity Flow Rate† (mL/hr) |
|-------------------|----------------------|----------------------------|
| Distal (16 Ga.) | 0.42 | 2656 |
| Medial (18 Ga.) | 0.38 | 1353 |
| Proximal (18 Ga.) | 0.41 | 1477 |

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

EU Authorized Representative and Importer:

Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL051217 R01 (2024-01)

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number



7 Fr. 3 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

REF BR-22703

Arrowgard Blue® Three-Lumen CVC

