

Dimensions:
Width: 6-9/64" (15.6 cm)
Length: 8-21/32" (22 cm)

PB-50014-100B
PBZ-50014-103A

T50014100A6 Dimensions: 5 7/8" x 7 5/8"



REF CV-20854
(IPNXXXXXX)



Arrowg+ard Blue® Four-Lumen CVC

Contents:

- 1: Arrowg+ard Blue® Four-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 11 cm
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 4: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F). Fluid path components are non-pyrogenic.
Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.

Language selection buttons: en, be, bg, hr, cs, da, nl, et, fi, fr, de, el, hu, is, it, lv, lt, no, pl, pt, ro, ru, sr, sk, sl, es, sv, tr, uk

MD (Medical Device) icon

STERILE icon

Latex-free icon (LATEX with slash)

Temperature icon (25°C / 77°F)



Lumen

Priming
Volume*
(mL)

Gravity
Flow
Rate†
(mL/hr)

Distal (16 Ga.)	0.37	3063
Medial 1 (14 Ga.)	0.43	5445
Medial 2 (18 Ga.)	0.35	838
Proximal (18 Ga.)	0.38	946

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

CAS # 7440-48-4

EU Authorized Representative and Importer:

Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Dublin Co. Westmeath, Ireland

CE 2797

All components are CE 2797 unless otherwise noted.

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

LBL051506 R03 (2024-09)



8.5 Fr. 4 Lumen

11 cm catheter length

.032 inch dia. spring-wire guide

REF CV-20854

Arrowg+ard Blue® Four-Lumen CVC

