

Dimensions:
 Width: 6-9/64" (15.6 cm)
 Length: 8-21/32" (22 cm)

PB-50014-100B
PBZ-50014-103A

T50014100A6 Dimensions: 5 7/8" x 7 5/8"



REF CV-24701-E
 (IPNXXXXXX)

ARROW

Arrowgard Blue® One-Lumen CVC

Contents:

- 1: Arrowgard Blue® One-Lumen Catheter: 14 Ga. (2.2 mm OD) x 20 cm Origin: XX XXXX XX XXX XXXXXXXX XXX XXXXXXXX XXXXXXXXXX
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Pressure Transduction Probe Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: SecondSite™ Adjustable Hub: Fastener Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX
- 1: SecondSite™ Adjustable Hub: Catheter Clamp Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXX


All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
 California Prop. 65

WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F). Fluid path components are non-pyrogenic.

Contraindications: The Arrowgard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (14 Ga.)	0.57	4905

* Priming volumes are approximate and are done without accessories.
 † Flow rate values are approximate and are determined using deionized water at 100 cm head height.

CAS # 7440-48-4
 EU Authorized Representative and Importer:
EC REP Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 2797 Co. Westmeath, Ireland
All components are CE 2797 unless otherwise noted.

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Set Product of xxx
 Packaged in xxx

Teleflex®
www.teleflex.com

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXXX
 (17)YYMMDD
 (11)YYMMDD
 (10)Lot Number

LBL051518 R03 (2025-09)



14 Ga. **1** Lumen

20 cm catheter length

.032 inch dia. spring-wire guide

REF CV-24701-E

Arrowgard Blue® One-Lumen CVC

ARROW