



REF ASK-14854-KR

8.5 Fr. 4 Lumen 30 cm catheter length .032 inch dia. spring-wire guide

# Four-Lumen CVC

Contents:

- 1: Four-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 30 cm

1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer

1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle

1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)

1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)

1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)

1: Introducer Needle: 18 Ga. x 1-1/2" (3.81 cm) XTW

1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW

1: Pressure Transduction Probe

1: Introducer Needle: 22 Ga. x 2-3/4" (7 cm) TW

1: Syringe: 10 mL Luer-Slip

1: Syringe: 5 mL Luer-Slip

1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL

1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm

4: Dust Cap: Non-Vented

1: SecondSite™ Adjustable Hub: Fastener

1: SecondSite™ Adjustable Hub: Catheter Clamp

2: Medicine Cup: 60 mL
- 1: SharpsAway® Disposal Cup

1: Drape: 35" x 50" (89 cm x 127 cm) with 4" (10 cm) fenestration, with adhesive

1: Drape: 30" x 44" (76 cm x 112 cm)

1: Needle Holder

1: Scalpel: #11

4: Stopcock: 4-way

2: Prep Sponge Swab

6: Gauze Pad: 4" x 4" (10 cm x 10 cm)

1: Dressing: OpSite®<sup>1</sup> 4-3/4" x 4" (12 cm x 10 cm)

1: Suture: 3-0 Silk with Curved Needle

2: Tray: Prep

<sup>1</sup>A registered trademark of T.J. Smith & Nephew, Limited.

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**  
**Not made with natural rubber latex.**  
Sterilized using ethylene oxide unless otherwise indicated in the contents list.  
Fluid path components are non-pyrogenic.



Lumen

Priming  
Volume\*  
(mL)

Gravity  
Flow  
Rate†  
(mL/hr)

Distal (16 Ga.)	0.51	1826
Medial 1 (14 Ga.)	0.62	3303
Medial 2 (18 Ga.)	0.41	830
Proximal (18 Ga.)	0.44	863

\* Priming volumes are approximate and are done without accessories.  
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



LOT



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Teleflex®  
www.teleflex.com

LBL053679 R00 (2021-05)



REF ASK-14854-KR

8.5 Fr. 4 Lumen 30 cm catheter length .032 inch dia. spring-wire guide

# Four-Lumen CVC