



REF

BR-14853-EK

(IPNXXXXXX)

8.5 Fr.

3 Lumen

30 cm catheter length

.032 inch dia. spring-wire guide

Three-Lumen CVC




Contents:









- 1: Three-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 30 cm
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Scalpel: #11
- 3: Stopcock: 4-way High-Flow
- All components are CE 2797 unless otherwise noted.


Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.





CAS # 7440-48-4

en

be

bg

hr

cs

da

nl

et

fi

fr

de

el

hu

is

it

lv

lt

no

pl

pt

ro

ru

sr

sk


sl

es

sv

tr

uk

 Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.59	2432
Medial (14 Ga.)	1.1	4914
Proximal (16 Ga.)	0.59	2083

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



EU Authorized Representative and Importer:

EC REP

Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA


Product of xxx
Packaged in xxx

LOT

LotNumber



YYYY-MM-DD



YYYY-MM-DD



(01)XXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)LotNumber

Teleflex®

www.teleflex.com

LBL054617 R02 (2024-08)



8.5 Fr.

3 Lumen

30 cm catheter length

.032 inch dia. spring-wire guide

REF

BR-14853-EK

Three-Lumen CVC

