

Product of xxx  
Packaged in xxx



REF CK-03000-B (IPNXXXXXX)

7 Fr. 2 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Two-Lumen CVC

Contents:

- 1: Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 16 cm
  - 1: Spring-Wire Guide: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other)
  - 1: Catheter: 18 Ga. x 1-9/16" (3.97 cm) Radiopaque over 20 Ga. RW Introducer Needle
  - 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
  - 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 3 mL Luer-Lock Syringe
  - 1: Syringe: 3 mL Luer-Lock
  - 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
  - 2: Dust Cap: Non-Vented
  - 1: SharpsAway® Disposal Cup
  - 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
  - 1: Scalpel: #11
  - 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
  - 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
  - 1: Suture: 3-0 Braided Silk with Straight Needle
- All components are CE 2797 unless otherwise noted.

Rx only  
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
Not made with natural rubber latex.  
Fluid path components are non-pyrogenic.

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.4	2681
Proximal (16 Ga.)	0.5	2396

\* Priming volumes are approximate and are done without accessories.  
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

en be bg hr cs da nl et fi fr de el hu is it lv lt no pl pt  
ro ru sr sk sl es sv tr uk



LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

EU Authorized Representative and Importer:  
Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland  
Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA  
Teleflex®  
[www.teleflex.com](http://www.teleflex.com)  
LBL054692 R01 (2024-02)



REF CK-03000-B

7 Fr. 2 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Two-Lumen CVC

