



REF EU-14403-N

4 Fr. 3 Lumen

13 cm catheter length

.018 inch dia. spring-wire guide

# Three-Lumen CVC

**Contents:**

- 1: Three-Lumen Catheter: 4 Fr. (1.5 mm OD) x 13 cm
- 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
- 1: Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SharpsAway® II Locking Disposal Cup
- 1: Safety Scalpel: #11

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)**  
**Not made with natural rubber latex.**  
 Fluid path components are non-pyrogenic.

All components are CE 2797 unless otherwise noted.



CAS # 7440-48-4



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (20 Ga.)	0.30	482
Medial (23 Ga.)	0.31	86
Proximal (23 Ga.)	0.34	99

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

EU Authorized Representative and Importer:

EC REP Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 2797 Co. Westmeath, Ireland

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL054725 R01 (2023-09)

LOT



4 Fr. 3 Lumen

13 cm catheter length

.018 inch dia. spring-wire guide

REF EU-14403-N

# Three-Lumen CVC