



REF EU-14553-N

(IPNXXXXXX)

5.5 Fr. 3 Lumen 30 cm catheter length .018 inch dia. spring-wire guide

Three-Lumen CVC for Femoral Vein Insertion

Contents:

- 1: Three-Lumen Catheter: 5.5 Fr. (1.9 mm OD) x 30 cm
- 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
- 1: Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway® II Locking Disposal Cup
- 1: Safety Scalpel: #11

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex.
 Fluid path components are non-pyrogenic.

All components are CE 2797 unless otherwise noted.



CAS # 7440-48-4

en be bg hr cs da nl et fi fr de el hu is it lv lt no
 pl pt ro ru sr sk sl es sv tr uk



Lumen

Priming Volume* (mL)

Gravity Flow Rate† (mL/hr)

Distal (20 Ga.)	0.39	244
Medial (22 Ga.)	0.38	211
Proximal (22 Ga.)	0.41	211

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



EU Authorized Representative and Importer:



Teleflex Medical

IDA Business and

Technology Park

Dublin Road, Athlone

Co. Westmeath, Ireland



All components are CE 2797 unless otherwise noted.

Arrow International LLC

Subsidiary of Teleflex Incorporated

3015 Carrington Mill Blvd.

Morrisville, NC 27560 USA

Product of xxx
 Packaged in xxx

LOT LotNumber

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXXX
 (17)YYMMDD
 (11)YYMMDD
 (10)LotNumber

Teleflex®
www.teleflex.com

LBL054807 R02 (2024-08)



5.5 Fr. 3 Lumen 30 cm catheter length .018 inch dia. spring-wire guide

REF EU-14553-N

Three-Lumen CVC for Femoral Vein Insertion

