



REF EU-17702-N

(IPNXXXXXX)

7 Fr. 2 Lumen

20 cm catheter length

.032 inch dia. spring-wire guide

Two-Lumen CVC

Contents:

- 1: Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm
  - 1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
  - 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
  - 1: Syringe: 5 mL Luer-Slip
  - 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
  - 1: SecondSite™ Adjustable Hub: Fastener
  - 1: SecondSite™ Adjustable Hub: Catheter Clamp
- All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
Not made with natural rubber latex.  
Fluid path components are non-pyrogenic.



CAS # 7440-48-4

en be bg hr cs da nl et fi fr de el hu is it lv lt no  
pl pt ro ru sr sk sl es sv tr uk



Lumen

Priming  
Volume\*  
(mL)

Gravity  
Flow  
Rate†  
(mL/hr)

Distal (14 Ga.)	0.60	5337
Proximal (18 Ga.)	0.43	1453

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



EU Authorized  
Representative and  
Importer:



Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland



All components are CE 2797 unless  
otherwise noted.



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Product of xxx  
Packaged in xxx

LOT LotNumber



YYYY-MM-DD



YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)LotNumber

Teleflex®  
[www.teleflex.com](http://www.teleflex.com)

LBL055066 R02 (2024-08)



7 Fr. 2 Lumen

20 cm catheter length

.032 inch dia. spring-wire guide

REF EU-17702-N

Two-Lumen CVC

