



REF ES-24402 (IPNXXXXXX)

4 Fr. 2 Lumen 13 cm catheter length .018 inch dia. spring-wire guide

Pediatric Arrowg+ard Blue® Two-Lumen CVC

Contents:

- 1: Arrowg+ard Blue® Two-Lumen Catheter: 4 Fr. (1.5 mm OD) x 13 cm
 - 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends)
 - 1: Introducer Needle: 22 Ga. x 1-1/16" (2.67 cm) TW
 - 1: Syringe: 3 mL Luer-Slip
 - 1: Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
 - 1: Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
 - 2: Dust Cap: Non-Vented
 - 1: SecondSite™ Adjustable Hub: Catheter Clamp
 - 1: SecondSite™ Adjustable Hub: Fastener
- All components are CE 2797 unless otherwise noted.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F). Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.



CAS # 7440-48-4



| Lumen | Priming Volume* (mL) | Gravity Flow Rate† (mL/hr) |
|-------------------|----------------------|----------------------------|
| Distal (22 Ga.) | 0.26 | 422 |
| Proximal (22 Ga.) | 0.30 | 191 |

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

EU Authorized Representative and Importer:

Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Product of xxx
 Packaged in xxx

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXXX
 (17)YYMMDD
 (11)YYMMDD
 (10)Lot Number

Teleflex®
 www.teleflex.com

LBL055201 R02 (2024-10)



4 Fr. 2 Lumen 13 cm catheter length .018 inch dia. spring-wire guide

REF ES-24402

Pediatric Arrowg+ard Blue® Two-Lumen CVC

