



REF AK-04150-E-S (IPNXXXXXX)



20 Ga. 1 Lumen 12 cm catheter length .025 inch dia. spring-wire guide

Small Vein Jugular Puncture

Contents:

- 1: One-Lumen Catheter: 20 Ga. (1.2 mm OD) x 12 cm Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 2: Syringe: 3 mL Luer-Lock Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Tissue Dilator: 4 Fr. (1.35 mm) x 7.6 cm Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Tissue Dilator: 4 Fr. (1.35 mm) x 4.1 cm Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: SecondSite™ Adjustable Hub: Catheter Clamp Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: SecondSite™ Adjustable Hub: Fastener Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Scalpel: #11 Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Prep Sponge Swab Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 3: Gauze Pad: 2" x 2" (5 cm x 5 cm) Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX

All components are CE 2797 unless otherwise noted.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



Lumen

Priming Volume* (mL)


Gravity Flow Rate† (mL/hr)

Distal (20 Ga.)	0.1	1531
-----------------	-----	------

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

- en be bg hr cs da nl et fi fr de el hu is it lv lt no pl pt
- ro ru sr sk sl es sv tr uk



EU Authorized Representative and Importer:



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Set Product of xxx
Packaged in xxx



www.teleflex.com
LBL055584 R03 (2024-11)



REF AK-04150-E-S

20 Ga. 1 Lumen

12 cm catheter length

.025 inch dia. spring-wire guide

Small Vein Jugular Puncture

