

Dimensions:  
Width: 6-9/64" (15.6 cm)  
Length: 8-21/32" (22 cm)

PB-50014-100B  
PBZ-50014-103A

T50014100A6 Dimensions: 5 7/8" x 7 5/8"



REF CV-12712  
(IPNXXXXXX)

ARROW®


# Two-Lumen CVC

## Contents:

- 1: Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 16 cm
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
Not made with natural rubber latex. Fluid path components are non-pyrogenic.



Lumen

Priming  
Volume\*  
(mL)

Gravity  
Flow  
Rate†  
(mL/hr)

Distal (14 Ga.)	0.52	5212
Proximal (18 Ga.)	0.38	1552

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



CAS # 7440-48-4



EU Authorized  
Representative and  
Importer:



Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
2797 Co. Westmeath, Ireland



All components are CE 2797  
unless otherwise noted.



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Product of xxx  
Packaged in xxx

Teleflex®  
[www.teleflex.com](http://www.teleflex.com)



Lot Number



YYYY-MM-DD



YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

LBL055599 R03 (2024-11)



7 Fr. 2 Lumen

16 cm  
catheter  
length

.032 inch dia.  
spring-wire  
guide

REF CV-12712

Two-Lumen CVC

ARROW®