



REF ES-04150

20<sup>Ga.</sup> 1<sup>Lumen</sup> 12<sup>cm</sup> catheter length .025<sup>inch dia.</sup> spring-wire guide

# Pediatric Jugular Puncture

**Contents:**

- 1: One-Lumen Catheter: 20 Ga. (1.2 mm OD) x 12 cm
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW
- 1: Syringe: 3 mL Luer-Slip
- 1: Tissue Dilator: 4 Fr. (1.35 mm) x 7.6 cm
- 1: Tissue Dilator: 4 Fr. (1.35 mm) x 4.1 cm
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener

All components are CE 2797 unless otherwise noted.

**Rx only**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU**

**Not made with natural rubber latex.**  
Fluid path components are non-pyrogenic.



CAS # 7440-48-4

- en
- be
- bg
- hr
- cs
- da
- nl
- et
- fi
- fr
- de
- el
- hu
- is
- it
- lv
- lt
- no
- pl
- pt
- ro
- ru
- sr
- sk
- sl
- es
- sv
- tr
- uk

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (20 Ga.)	0.1	1531

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

EU Authorized Representative and Importer

Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
www.teleflex.com

LBL055763 R01 (2023-09)

LOT



20<sup>Ga.</sup> 1<sup>Lumen</sup> 12<sup>cm</sup> catheter length .025<sup>inch dia.</sup> spring-wire guide

REF ES-04150

# Pediatric Jugular Puncture