

Product of XXX  
Packaged in XXX



REF **ES-14402** (IPNXXXXXX)

**4** Fr. **2** Lumen **13** cm catheter length **.018** inch dia. spring-wire guide

# Pediatric Two-Lumen CVC

**Contents:**

- 1: Two-Lumen Catheter: 4 Fr. (1.5 mm OD) x 13 cm
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends)
- 1: Introducer Needle: 22 Ga. x 1-1/16" (2.67 cm) TW
- 1: Syringe: 3 mL Luer-Slip
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener

**Rx only**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)**  
**California Prop. 65**  
 ⚠ **WARNING:** Cancer and Reproductive Harm.  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)  
**Not made with natural rubber latex.**  
 Fluid path components are non-pyrogenic.

**All components are CE 2797 unless otherwise noted.**



CAS # 7440-48-4



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (22 Ga.)	0.27	363
Proximal (22 Ga.)	0.31	277

\* Priming volumes are approximate and are done without accessories.  
 † Flow rate values are approximate and are determined using deionized water at 100 cm head height.

EU Authorized Representative and Importer:

**EC REP** Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland  
**CE** 2797

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL055793 R01 (2024-01)

**LOT** LotNumber

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXX  
 (17)YYMMDD  
 (11)YYMMDD  
 (10)LotNumber



**4** Fr. **2** Lumen **13** cm catheter length **.018** inch dia. spring-wire guide

REF ES-14402

# Pediatric Two-Lumen CVC

