



REF CV-12702

## Two-Lumen CVC

**Contents:**



- 1: Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 16 cm
- 1: Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 2797 unless otherwise noted.

**Rx only Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU) California Prop. 65

**⚠ WARNING:** Cancer and Reproductive Harm. [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**Not made with natural rubber latex.** Fluid path components are non-pyrogenic.



| Lumen             | Priming Volume* (mL) | Gravity Flow Rate† (mL/hr) |
|-------------------|----------------------|----------------------------|
| Distal (18 Ga.)   | 0.34                 | 1220                       |
| Proximal (14 Ga.) | 0.51                 | 4548                       |

\* Priming volumes are approximate and are done without accessories.  
 † Flow rate values are approximate and are determined using deionized water at 100 cm head height.

CAS # 7440-48-4

EU Authorized Representative and Importer

EC REP Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone Co. Westmeath, Ireland

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL055814 R01 (2023-09)

LOT



7 Fr. 2 Lumen 16 cm catheter length .025 inch dia. spring-wire guide

REF CV-12702

## Two-Lumen CVC