



REF CS-16553-E

(IPNXXXXXX)

5.5 Fr. 3 Lumen

13 cm catheter length

.018 inch dia. spring-wire guide

Pediatric Three-Lumen CVC

Contents:

- 1: Three-Lumen Catheter: 5.5 Fr. (1.9 mm OD) x 13 cm
  - 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
  - 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
  - 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
  - 1: Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
  - 1: Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
  - 3: Dust Cap: Non-Vented
  - 1: SecondSite™ Adjustable Hub: Fastener
  - 1: SecondSite™ Adjustable Hub: Catheter Clamp
- All components are CE 2797 unless otherwise noted.

Rx only  
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
Not made with natural rubber latex.  
Fluid path components are non-pyrogenic.

CAS # 7440-48-4

en

be

bg

hr

cs

da

nl

et

fi

fr

de

el

hu

is

it

lv

lt

no

pl

pt

ro

ru

sr

sk

sl

es

sv

tr

uk

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (20 Ga.)	0.32	988
Medial (22 Ga.)	0.38	370
Proximal (22 Ga.)	0.34	422

\* Priming volumes are approximate and are done without accessories.  
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

EU Authorized Representative and Importer:

Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Product of xxx  
Packaged in xxx

LOT LotNumber

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)LotNumber

Teleflex®  
[www.teleflex.com](http://www.teleflex.com)

LBL056000 R02 (2024-05)



5.5 Fr. 3 Lumen

13 cm catheter length

.018 inch dia. spring-wire guide

REF CS-16553-E

Pediatric Three-Lumen CVC

