

Dimensions:
 Width: 6-9/64" (15.6 cm)
 Length: 8-21/32" (22 cm)

PB-50014-100B
PBZ-50014-103A

T50014100A6 Dimensions: 5 7/8" x 7 5/8"



REF CV-15553
 (IPNXXXXXX)

ARROW


Pediatric Three-Lumen CVC

Contents:

- 1: Three-Lumen Catheter: 5.5 Fr. (1.9 mm OD) x 8 cm
- 1: Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
 Not made with natural rubber latex. Fluid path components are non-pyrogenic.



Lumen

Priming Volume* (mL)

Gravity Flow Rate† (mL/hr)

Distal (20 Ga.)	0.29	1253
Medial (22 Ga.)	0.31	601
Proximal (22 Ga.)	0.34	872

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



CAS # 7440-48-4



EU Authorized Representative and Importer:



Teleflex Medical
 IDA Business and
 Technology Park



Dublin Road, Athlone
 2797 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.



Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Product of xxx
 Packaged in xxx

Teleflex
 www.teleflex.com



Lot Number



YYYY-MM-DD



YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXXX
 (17)YYMMDD
 (11)YYMMDD
 (10)Lot Number

LBL056003 R03 (2024-10)



5.5 Fr. **3** Lumen

8 cm catheter length

.018 inch dia. spring-wire guide

REF CV-15553

Pediatric Three-Lumen CVC

ARROW