



**REF CV-15553**

## Pediatric Three-Lumen CVC



**Contents:**


- 1: Three-Lumen Catheter: 5.5 Fr. (1.9 mm OD) x 8 cm
  - 1: Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
  - 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
  - 1: Syringe: 5 mL Luer-Slip
  - 1: Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
  - 1: Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
  - 3: Dust Cap: Non-Vented
  - 1: SecondSite™ Adjustable Hub: Fastener
  - 1: SecondSite™ Adjustable Hub: Catheter Clamp
- All components are CE 2797 unless otherwise noted.**

**Rx only Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU) California Prop. 65


**⚠ WARNING:** Cancer and Reproductive Harm. [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

**Not made with natural rubber latex.** Fluid path components are non-pyrogenic.



 Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (20 Ga.)	0.29	1253
Medial (22 Ga.)	0.31	601
Proximal (22 Ga.)	0.34	872

\* Priming volumes are approximate and are done without accessories.  
 † Flow rate values are approximate and are determined using deionized water at 100 cm head height.

 CAS # 7440-48-4

 EU Authorized Representative and Importer:

**EC REP** Teleflex Medical  
 IDA Business and  
 Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland  
**CE** 2797

**Arrow International LLC**  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL056004 R01 (2023-10)

**LOT**



**5.5** Fr. **3** Lumen

**8** cm catheter length

**.018** inch dia. spring-wire guide

**REF CV-15553**

**Pediatric Three-Lumen CVC**