

Product of XXX  
Packaged in XXX



REF ES-04218 (IPNXXXXXX)

18<sup>Ga.</sup> 1<sup>Lumen</sup> 20<sup>cm</sup> catheter length .025<sup>inch dia.</sup> spring-wire guide

# One-Lumen CVC

## Contents:

- 1: One-Lumen Catheter: 18 Ga. (1.5 mm OD) x 20 cm
- 1: Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener

## Rx only

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**  
[www.teleflex.com/IFU](http://www.teleflex.com/IFU)

**Not made with natural rubber latex.**

Fluid path components are non-pyrogenic.

All components are CE 2797 unless otherwise noted.



CAS # 7440-48-4

- en be bg hr cs da nl et fi fr
- de el hu is it lv lt no pl pt
- ro ru sr sk sl es sv tr uk



Lumen

Priming Volume\* (mL)

Gravity Flow Rate† (mL/hr)

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (18 Ga.)	0.21	2025

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



EU Authorized Representative and Importer:



Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA



LBL056093 R01 (2024-01)

LOT LotNumber

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)LotNumber



18<sup>Ga.</sup> 1<sup>Lumen</sup> 20<sup>cm</sup> catheter length .025<sup>inch dia.</sup> spring-wire guide



REF ES-04218

# One-Lumen CVC