



REF **CS-25703-E** (IPNXXXXXX)

**7** Fr. **3** Lumen **20** cm catheter length **.032** inch dia. spring-wire guide

# Arrowg+ard Blue® Three-Lumen CVC

**Contents:**

- 1: Arrowg+ard Blue® Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm
  - 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
  - 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
  - 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
  - 1: Pressure Transduction Probe
  - 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
  - 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
  - 3: Dust Cap: Non-Vented
  - 1: SecondSite™ Adjustable Hub: Fastener
  - 1: SecondSite™ Adjustable Hub: Catheter Clamp
- All components are CE 2797 unless otherwise noted.**

**Rx only**

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)

Not made with natural rubber latex.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F). Fluid path components are non-pyrogenic.

**Contraindications:** The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.



Lumen


Priming Volume\* (mL)

Gravity Flow Rate† (mL/hr)

Distal (16 Ga.)	0.4	2258
Medial (18 Ga.)	0.4	1046
Proximal (18 Ga.)	0.4	1154

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



EU Authorized Representative and Importer:  
**EC REP** Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland  
**CE** 2797

**All components are CE 2797 unless otherwise noted.**

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**LOT** LotNumber

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXXX  
 (17)YYMMDD  
 (11)YYMMDD  
 (10)LotNumber

Product of xxx  
 Packaged in xxx

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL056423 R02 (2024-07)



REF **CS-25703-E**

**7** Fr. **3** Lumen

**20** cm catheter length

**.032** inch dia. spring-wire guide



# Arrowg+ard Blue® Three-Lumen CVC