

Dimensions:  
 Width: 8-3/4" (22,2 cm)  
 Length: 11-1/8" (28,3 cm)

**PB-10001-000A**  
**PBZ-10001-000A**

T-12703-100C  
 Width: 8.7598 in = 222.5mm  
 Length: 10.374 in = 263.5mm

**ARROW**



**REF** EU-15703-CVTS (IPNXXXXXX)

**7** Fr. **3** Lumen **20** cm catheter length **.032** inch dia. spring-wire guide

# Three-Lumen CVC

**Contents:**

- 1: Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 1: ECG Extension Cable
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway® II Locking Disposal Cup
- 1: Safety Scalpel: #11

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**  
[www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
**Not made with natural rubber latex.**  
 Fluid path components are non-pyrogenic.

All components are CE 2797 unless otherwise noted.

CAS # 7440-48-4

en be bg hr cs da nl et fi fr de el hu is it lv lt no pl pt ro ru sr  
 sk sl es sv tr uk

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.47	2656
Medial (18 Ga.)	0.40	1245
Proximal (18 Ga.)	0.42	1328

\* Priming volumes are approximate and are done without accessories.  
 † Flow rate values are approximate and are determined using deionized water at 100 cm head height.

EU Authorized Representative and Importer:  
**Teleflex Medical**  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland  
**CE 2797**  
 All components are CE 2797 unless otherwise noted.

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**LOT** LotNumber  
 YYYY-MM-DD YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
 (17)YYMMDD  
 (11)YYMMDD  
 (10)LotNumber

Product of xxx  
 Packaged in xxx

**Teleflex**  
[www.teleflex.com](http://www.teleflex.com)  
 LBL056836 R02 (2024-08)



# Three-Lumen CVC

**REF** EU-15703-CVTS

**7** Fr. **3** Lumen **20** cm catheter length **.032** inch dia. spring-wire guide

**ARROW**