

Dimensions:  
 Width: 6-9/64" (15.6 cm)  
 Length: 8-21/32" (22 cm)

PB-50014-100B  
**PBZ-50014-103A**

T50014100A6 Dimensions: 5 7/8" x 7 5/8"



**REF** CV-50014  
 (IPNXXXXXX)

**ARROW**


# One-Lumen CVC

**Contents:**

- 1: One-Lumen Catheter: 14 Ga. (2.2 mm OD) x 20 cm, Side Holes
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 2797 unless otherwise noted.

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU**  
 Not made with natural rubber latex. Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (14 Ga.)	0.4	6748

\* Priming volumes are approximate and are done without accessories.  
 † Flow rate values are approximate and are determined using deionized water at 100 cm head height.

CAS # 7440-48-4

EU Authorized Representative and Importer:

**EC REP** Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 2797 Co. Westmeath, Ireland  
**CE** All components are CE 2797 unless otherwise noted.

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

Product of xxx  
 Packaged in xxx

**Teleflex**  
 www.teleflex.com

**LOT** Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXX  
 (17)YYMMDD  
 (11)YYMMDD  
 (10)Lot Number

LBL056954 R03 (2024-12)



**14** Ga. **1** Lumen

**20** cm catheter length

**.032** inch dia. spring-wire guide

**REF** CV-50014

**One-Lumen CVC**

**ARROW**