



REF

EU-15955-CVT

(IPNXXXXXX)

9.5Fr.

5Lumen

20cm catheter length

.032inch dia. spring-wire guide

Five-Lumen CVC

Contents:

- 1: Five-Lumen Catheter: 9.5 Fr. (3.3 mm OD) x 20 cm
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
- 1: ECG Extension Cable
- 5: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Scalpel: #11

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



CAS # 7440-48-4

enbebghrcsdaneletfifrdeelhuisisitlvltno

plptroaruskslsesvstruk



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.48	2988
Medial 1 (14 Ga.)	0.57	4474
Medial 2 (18 Ga.)	0.38	1038
Medial 3 (18 Ga.)	0.38	963
Proximal (18 Ga.)	0.41	1204

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



EU Authorized Representative and Importer:



Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland



All components are CE 2797 unless otherwise noted.



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Product of xxx  
Packaged in xxx

LOT Lot Number



YYYY-MM-DD



YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

Teleflex®  
[www.teleflex.com](http://www.teleflex.com)

LBL057152 R02 (2025-03)



9.5Fr.

5Lumen

20cm catheter length

.032inch dia. spring-wire guide

REF EU-15955-CVT

Five-Lumen CVC

