



REF ML-00701

One-Lumen CVC for use only with Arrow MAC™ Two-Lumen Central Venous Access Device


Contents:

1: One-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm

1: Dust Cap: Non-Vented

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU Not made with natural rubber latex. Fluid path components are non-pyrogenic.



STERILE EO



Lumen

Priming
Volume*
(mL)

Gravity
Flow
Rate†
(mL/hr)

Distal (12 Ga.)

0.66

7727

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



CAS # 7440-48-4



EU Authorized
Representative and
Importer:

EC REP

CE
2797

Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL057176 R02 (2023-10)

LOT



7 Fr.

1 Lumen

14 cm catheter
length with
MAC

REF ML-00701

One-Lumen CVC for use only with Arrow MAC™
Two-Lumen Central Venous Access Device