

Dimensions:
Width: 6-9/64" (15,6 cm)
Length: 8-21/32" (22 cm)

PB-50014-100B
PBZ-50014-103A

T50014100A6 Dimensions: 5 7/8" x 7 5/8"



REF ML-00701
(IPNXXXXXX)

ARROW®


One-Lumen CVC for use only with Arrow MAC™ Two-Lumen
Central Venous Access Device

Contents:

- 1: One-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm
1: Dust Cap: Non-Vented

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use.
Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex. Fluid path components are non-pyrogenic.



Lumen

Priming
Volume*
(mL)

Gravity
Flow
Rate†
(mL/hr)

Distal (12 Ga.) 0.66 7727

* Priming volumes are approximate and are done without accessories.
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

CAS # 7440-48-4

EU Authorized
Representative and
Importer:

EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
2797 Co. Westmeath, Ireland
All components are CE 2797
unless otherwise noted.

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com

LOT Lot Number



YYYY-MM-DD



YYYY-MM-DD



(01)XXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

LBL057176 R03 (2024-09)



7 Fr. 1 Lumen 14 cm catheter
length with MAC

REF ML-00701

One-Lumen CVC for use only with Arrow MAC™
Two-Lumen Central Venous Access Device

ARROW®