



REF CV-12122-FBZ

12 Fr. 2 Lumen 16 cm catheter length .035 inch dia. spring-wire guide

# Two-Lumen Hemodialysis Catheter for High Volume Infusions

**Contents:**

- 1: Two-Lumen Catheter: 12 Fr. (4.06 mm OD) x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 2: Dust Cap: Non-Vented

All components are CE 2797 unless otherwise noted.

**Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU California Prop. 65**

**⚠ WARNING:** Cancer and Reproductive Harm. www.P65Warnings.ca.gov

**Not made with natural rubber latex.**

Fluid path components are non-pyrogenic.

**Contraindications: The large-bore two-lumen catheter is not designed for long-term (≥ 30 days) hemodialysis or for use in patients with thrombosed vessels.**



| Lumen             | Priming Volume* (mL) | Flow Rate (mL/min) | Venous Pressure Setting†† (mmHg) |
|-------------------|----------------------|--------------------|----------------------------------|
| Distal (12 Ga.)   | 1.4                  | 220                | 100                              |
|                   |                      | 370                | 200                              |
| Proximal (12 Ga.) | 1.3                  | 220                | 100                              |
|                   |                      | 370                | 200                              |

\* Priming volumes are approximate and are done without accessories.  
 †† Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an in vitro experiment using blood analog solution.

EU Authorized Representative and Importer:  
 Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 2797 Co. Westmeath, Ireland

Arrow International LLC  
 Subsidiary of Teleflex Incorporated  
 3015 Carrington Mill Blvd.  
 Morrisville, NC 27560 USA

**Teleflex®**  
 www.teleflex.com

LBL057798 R00 (2021-05)

LOT

