



REF CS-12142-F (IPNXXXXXX)

14^{Fr.} 2^{Lumen} 15^{cm catheter length} .038^{inch dia. spring-wire guide}


Two-Lumen Hemodialysis Catheter for High Volume Infusions

Contents:

- 1: Two-Lumen Catheter: 14 Fr. (4.72 mm OD) x 15 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .038" (0.97 mm) dia. x 27-1/2" (70 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 14 Fr. (4.7 mm) x 14 cm
- 1: Tissue Dilator 16 Fr. (5.3 mm) x 4 1/2" (11.4 cm)
- 2: Dust Cap: Non-Vented
- 1: Drape: 24" x 36" (61 cm x 91 cm)
- 1: Scalpel: #11

All components are CE 2797 unless otherwise noted.

Rx only
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex.
Fluid path components are non-pyrogenic.
Contraindications: The large-bore two-lumen catheter is not designed for long-term (≥ 30 days) hemodialysis or for use in patients with thrombosed vessels.

|  Lumen | Priming Volume* (mL) | Flow Rate (mL/min) | Venous Pressure Setting†† (mmHg) |
|---|----------------------|--------------------|----------------------------------|
| Distal (10 Ga.) | 1.4 | 250 400 | 70 100 |
| Proximal (10 Ga.) | 1.2 | 250 400 | 70 150 |

* Priming volumes are approximate and are done without accessories.
†† Venous Pressure settings are approximate and were determined by varying the blood flow rates from 250 to 400 mL/min in an in vitro experiment using blood analog solution.

- enbgcsdaneletfifrdeelhuitletonoplastro
- skslsvtr





LOT Lot Number
YYYY-MM-DD YYYY-MM-DD



(01)XXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Product of xxx
Packaged in xxx

EU Authorized Representative and Importer:  Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
All components are CE 2797 unless otherwise noted.

 Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA


www.teleflex.com
LBL057824 R04 (2024-11)



Two-Lumen Hemodialysis Catheter for High Volume Infusions



REF CS-12142-F 14^{Fr.} 2^{Lumen} 15^{cm catheter length} .038^{inch dia. spring-wire guide}