

Product of xxx
Packaged in xxx



REF **AD-09903-S** (IPNXXXXXX)

9 Fr. **10** cm sheath length

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 8 Fr. Catheters

Contents:

- 1: Arrow-Flex® Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
- 1: Dust Cap: Non-Vented
- 1: Stopcock: 4-way High-Flow

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

en fr de it pt es sv



LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

EU Authorized Representative and Importer:
EC REP Teleflex Medical
 IDA Business and
 Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Teleflex
 www.teleflex.com

LBL057930 R01 (2024-01)



REF **AD-09903-S**

9 Fr. **10** cm sheath length

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 8 Fr. Catheters

