

DO NOT BEND

ARROW®



REF SAC-00324

(IPNXXXXXX)

24^{Ga.} 2.5^{cm catheter length}

Arterial Catheterization Set

Contents:


- 1: Catheter: 24 Ga. x 2.5 cm Radiopaque Polyurethane
 1: Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 9-3/4" (25 cm) (Straight Soft Tip on Both Ends)
 1: Catheter: 24 Ga. x 3/4" (1.9 cm) Radiopaque over 26 Ga. Introducer Needle
All components are CE 2797 unless otherwise noted.

Catheter Priming Volume: 0.14 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.
www.teleflex.com/IFU

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



EU Authorized Representative and Importer:

EC REP



Teleflex Medical
 IDA Business and
 Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD

Product of xxx
Packaged in xxx

(01)XXXXXXXXXXXXXX
 (17)YYMMDD
 (11)YYMMDD
 (10)Lot Number

Teleflex®
www.teleflex.com

LBL058176 R02 (2025-03)



REF SAC-00324

24^{Ga.} 2.5^{cm catheter length}

Arterial Catheterization Set

ARROW®

DO NOT BEND