

DO NOT BEND



REF SAC-00324-PBX

(IPNXXXXXX)

24 Ga. 2.5 cm catheter length

Arterial Catheterization Set

Contents:

- 1: Catheter: 24 Ga. x 2.5 cm Radiopaque PEBA
 - 1: Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 9-3/4" (25 cm) (Straight Soft Tip on Both Ends)
 - 1: Catheter: 24 Ga. x 3/4" (1.9 cm) Radiopaque over 26 Ga. Introducer Needle
- All components are CE 2797 unless otherwise noted.**

Catheter Priming Volume: 0.14 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

www.teleflex.com/IFU

California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



- en bg cs da nl et fi fr de el hu
- it lv lt no pl pt ro ru sk sl es
- sv tr

Product of xxx
Packaged in xxx

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number



EU Authorized Representative and Importer:



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex
www.teleflex.com

LBL058179 R01 (2024-01)



REF SAC-00324-PBX

24 Ga. 2.5 cm catheter length

Arterial Catheterization Set



DO NOT BEND