

DO NOT BEND

ARROW®



REF SAC-00522

(IPNXXXXXX)

22^{Ga.}5^{cm}
catheter
length

Arterial Catheterization Set

Contents:

- 1: Catheter: 22 Ga. x 5 cm Radiopaque Polyurethane
 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
 1: Introducer Needle: 22 Ga. x 1-1/2" (4 cm) XTW


All components are CE 2797 unless otherwise noted.

Catheter Priming Volume: 0.16 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.
www.teleflex.com/IFU

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

EU Authorized
Representative and
Importer:

EC REP



Teleflex Medical
 IDA Business and
 Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD

Product of xxx
Packaged in xxx

(01)XXXXXXXXXXXXXX
 (17)YYMMDD
 (11)YYMMDD
 (10)Lot Number

Teleflex®
www.teleflex.com

LBL058188 R02 (2025-03)



REF SAC-00522

22^{Ga.}5^{cm}
catheter
length

Arterial Catheterization Set

ARROW®

DO NOT BEND