

DO NOT BEND



REF SAC-00522

(IPNXXXXXX)

22 Ga.

5 cm catheter length

# Arterial Catheterization Set

**Contents:**

- 1: Catheter: 22 Ga. x 5 cm Radiopaque Polyurethane
- 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
- 1: Introducer Needle: 22 Ga. x 1-1/2" (4 cm) XTW

All components are CE 2797 unless otherwise noted.

Catheter Priming Volume: 0.16 mL Priming volume is approximate.

**Rx only Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

[www.teleflex.com/IFU](http://www.teleflex.com/IFU)

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



- en bg cs da nl et fi fr de el hu
- it lv lt no pl pt ro ru sk sl es
- sv tr

Product of xxx  
Packaged in xxx

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

EU Authorized Representative and Importer:

EC REP Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Teleflex®  
[www.teleflex.com](http://www.teleflex.com)

LBL058188 R01 (2024-01)



REF SAC-00522

22 Ga.

5 cm catheter length

# Arterial Catheterization Set



DO NOT BEND