

DO NOT BEND



REF SAC-00524

(IPNXXXXXX)

24 Ga.

5 cm catheter length

Arterial Catheterization Set

Contents:

- 1: Catheter: 24 Ga. x 5 cm Radiopaque Polyurethane
- 1: Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 9-3/4" (25 cm) (Straight Soft Tip on Both Ends)
- 1: Catheter: 24 Ga. x 3/4" (1.9 cm) Radiopaque over 26 Ga. Introducer Needle

All components are CE 2797 unless otherwise noted.


Catheter Priming Volume: 0.15 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

www.teleflex.com/IFU

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



EU Authorized Representative and Importer:

EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD

Product of xxx
Packaged in xxx



(01)XXXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

Teleflex®
www.teleflex.com

LBL058194 R02 (2025-03)



REF SAC-00524

24 Ga.

5 cm catheter length

Arterial Catheterization Set



DO NOT BEND