

DO NOT BEND



REF SAC-01618

(IPNXXXXXX)

18 Ga. 16 cm catheter length

Arterial Catheterization Set

Contents:

- 1: Catheter: 18 Ga. x 16 cm Radiopaque Polyurethane
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW

All components are CE 2797 unless otherwise noted.

Catheter Priming Volume: 0.24 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

www.teleflex.com/IFU

California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



- en bg cs da nl et fi fr de el hu
- it lv lt no pl pt ro ru sk sl es
- sv tr

Product of xxx
Packaged in xxx

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

EU Authorized Representative and Importer:

EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL058236 R01 (2024-01)



REF SAC-01618

18 Ga. 16 cm catheter length

Arterial Catheterization Set



DO NOT BEND