

DO NOT BEND



REF SAC-01620

(IPNXXXXXX)

20 Ga. 16 cm catheter length

# Arterial Catheterization Set

**Contents:**

- 1: Catheter: 20 Ga. x 16 cm Radiopaque Polyurethane
- 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 19-3/4" (50 cm) (Straight Soft Tip on Both Ends) with Arrow Advancer
- 1: Introducer Needle: 20 Ga. x 2-3/4" (7 cm) XTW

All components are CE 2797 unless otherwise noted.


Catheter Priming Volume: 0.21 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

www.teleflex.com/IFU  
California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.  
www.P65Warnings.ca.gov

Not made with natural rubber latex.  
Fluid path components are non-pyrogenic.



EU Authorized Representative and Importer:

EC REP Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland  
2797

All components are CE 2797 unless otherwise noted.

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD

Product of xxx  
Packaged in xxx



(01)XXXXXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

Teleflex®  
www.teleflex.com

LBL058242 R02 (2025-03)



REF SAC-01620

20 Ga. 16 cm catheter length

# Arterial Catheterization Set



DO NOT BEND