



REF SI-09880-AG

8.5^{Fr.} 10^{cm} sheath length .035^{inch dia.} spring-wire guide

Arrowgard Blue® Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

Contents:

- 1: Arrowgard Blue® Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
 - 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
 - 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
 - 1: Obturator: 8 Fr.
 - 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
 - 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- All components are CE 2797 unless otherwise noted.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

Not made with natural rubber latex.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The Arrowgard Blue antimicrobial sheath introducer is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).

en cs fr de el it pl pt ru sl es sv tr



CAS # 7440-48-4



EU Authorized Representative and Importer:
 Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL058274 R01 (2023-09)