



REF IC-09880



8.5 Fr.

10 cm sheath length

.035 inch dia. spring-wire guide

PSI Set for use with 7 - 7.5 Fr. Catheters

Contents:

- 1: Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
- 1: Dust Cap: Non-Vented
- 1: Safety Scalpel: #11
- 1: Stopcock: 3-way High-Flow
- 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Caution: This Product Contains Dry Natural Rubber.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



LOT



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL058313 R00 (2021-11)



REF IC-09880

8.5 Fr.

10 cm sheath length

.035 inch dia. spring-wire guide

PSI Set for use with 7 - 7.5 Fr. Catheters

