



REF MTO-09801-CHA

8.5<sup>Fr.</sup>

10<sup>cm</sup>  
sheath  
length

.035<sup>inch dia.</sup>  
spring-wire  
guide



PSI Kit for use with 7 - 7.5 Fr. Catheters

Contents:

- 1: Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator

1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer

1: Catheter: 18 Ga. x 1-9/16" (3.97 cm) Radiopaque over 20 Ga. RW Introducer Needle

1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)

1: Injection Needle: 25 Ga. x 1" (2.54 cm)

1: Pressure Transduction Probe

1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe

1: Syringe: 3 mL Luer-Lock

1: Syringe: 5 mL Luer-Slip

1: Obturator: 8 Fr.

1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape

1: Dust Cap: Non-Vented

1: SharpsAway® Disposal Cup

1: Drape: 24" x 36" (61 cm x 91 cm)

1: Drape: 24" x 36" with 4" fenestration

1: Scalpel: #11

1: Stopcock: 3-way High-Flow
- 2: Prep Sponge Swab

2: Gauze Pad: 2" x 2" (5 cm x 5 cm)

5: Gauze Pad: 4" x 4" (10 cm x 10 cm)

1: Suture: 3-0 Nylon with Straight Needle

1: Tray: Prep

Rx only  
Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.  
Caution: This Product Contains Dry Natural Rubber.  
Sterilized using ethylene oxide unless otherwise indicated in the contents list.  
Fluid path components are non-pyrogenic.



Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA  
**Teleflex®**  
www.teleflex.com  
LBL058336 R01 (2022-02)



REF MTO-09801-CHA

8.5<sup>Fr.</sup>

10<sup>cm</sup>  
sheath  
length

.035<sup>inch dia.</sup>  
spring-wire  
guide

PSI Kit for use with 7 - 7.5 Fr. Catheters

