



REF ASK-09907-SB (IPNXXXXXX)

9 Fr.

10 cm sheath length

.025 inch dia. spring-wire guide

PSI Kit for use with 7 Fr. Catheters

Contents:

- 1: Arrow-Flex® Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator
 - 1: Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
 - 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
 - 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
 - 1: Pressure Transduction Probe
 - 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
 - 1: Injection Needle: SafetyGlide™¹ 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
 - 1: Injection Needle: SafetyGlide™¹ 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
 - 1: Obturator: 8 Fr.
 - 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
 - 1: Dust Cap: Non-Vented
 - 1: SharpsAway® II Locking Disposal Cup
 - 1: SharpsAway® Disposal Cup
 - 1: Maximal Barrier Drape™ with 4" fenestration
 - 1: Drape: 26" x 18" (66 cm x 46 cm) with adhesive
 - 1: Safety Scalpel: #11
 - 1: Stopcock: 4-way High-Flow
 - 2: Prep Sponge Swab
 - 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
 - 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
 - 1: Suture: 3-0 Braided Silk with Straight Needle
 - 1: Tray: Prep
- ¹A trademark of Becton, Dickinson and Company.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



LOT LotNumber

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)LotNumber

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of XXX
Packaged in XXX



LBL058350 R03 (2024-05)



REF ASK-09907-SB

9 Fr.

10 cm sheath length

.025 inch dia. spring-wire guide



PSI Kit for use with 7 Fr. Catheters

