

REF **ASK-09803-FM** (IPNXXXXXX)**8.5** Fr.**10** cm
sheath
length**.035** inch dia.
spring-wire
guide

PSI Kit for use with 7 - 7.5 Fr. Catheters



Contents:

- 1: Single-Lumen Infusion Catheter, SLIC®: 7 Fr. x 16 cm
 - 1: Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
 - 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
 - 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
 - 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
 - 1: Injection Needle: Eclipse™¹ 23 Ga. x 1-1/2" (3.81 cm)
 - 1: Injection Needle: Eclipse™¹ 25 Ga. x 1" (2.54 cm)
 - 1: Syringe: 3 mL Luer-Lock
 - 2: Syringe: 5 mL Luer-Slip
 - 1: Obturator: 8 Fr.
 - 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
 - 1: SharpsAway® Disposal Cup
 - 1: Maximal Barrier Drape™ with 4" fenestration
 - 1: Safety Scalpel: #11
 - 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
 - 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
 - 1: Suture: 3-0 Braided Silk with Straight Needle
- ¹A trademark of Becton, Dickinson and Company.

Not made with natural rubber latex.

