



REF ASK-00520-AHS (IPNXXXXXX)

20 Ga.

5 cm catheter length

.021 inch dia. spring-wire guide

Arterial Catheterization Kit

Contents:

- 1: Catheter: 20 Ga. x 5 cm Radiopaque PEBA
 - 1: Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
 - 1: Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
 - 1: Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
 - 1: Syringe: 5 mL Luer-Slip
 - 1: Dust Cap: Non-Vented
 - 1: SharpsAway® II Locking Disposal Cup
 - 1: Drape: 30" x 36" (76 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive
 - 1: Needle Holder
 - 1: Scissors
 - 3: Gauze Pad: 2" x 2" (5 cm x 5 cm)
 - 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
 - 1: Dressing: Tegaderm®² CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm)
 - 1: Suture: 3-0 Silk with Curved Needle
- ¹A trademark of Becton, Dickinson and Company.
²A registered trademark of 3M Company.

Catheter Priming Volume: 0.17 mL Priming volume is approximate.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



LOT LotNumber

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)LotNumber

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of XXX
Packaged in XXX



LBL059238 R02 (2025-03)



REF ASK-00520-AHS

20 Ga.

5 cm catheter length

.021 inch dia. spring-wire guide

Arterial Catheterization Kit

