



REF ASK-04220-HHC

20^{Ga.} 3.81^{cm} catheter length .018^{inch dia.} spring-wire guide

QuickFlash® Radial Artery Catheterization Kit

Contents:

- 1: QuickFlash® Radial Artery Catheter: 20 Ga. x 1-1/2" (3.81 cm) Radiopaque Polyurethane over 21 Ga. TW Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Wing Clip: Suture
- 1: SharpsAway® II Locking Disposal Cup
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: I.V. Advanced 6.5 cm x 7 cm

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



LOT



Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

LBL059286 R01 (2022-02)



REF ASK-04220-HHC

20^{Ga.} 3.81^{cm} catheter length .018^{inch dia.} spring-wire guide

QuickFlash® Radial Artery Catheterization Kit