



REF AK-04510-S (IPNXXXXXX)

20 Ga. 12 cm catheter length .025 inch dia. spring-wire guide

Arterial Line

Contents:

- 1: Indwelling Catheter: 20 Ga. x 12 cm Radiopaque FEP origin:xxxxxxxxxxxxxxxxxxxx
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW origin:xxxxxxxxxxxxxxxxxxxx
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 2: Syringe: 3 mL Luer-Lock origin:xxxxxxxxxxxxxxxxxxxx
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep origin:xxxxxxxxxxxxxxxxxxxx
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Solution with Hi-Lite Orange™ Tint origin:xxxxxxxxxxxxxxxxxxxx
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration origin:xxxxxxxxxxxxxxxxxxxx
- 1: Scalpel: #11 origin:xxxxxxxxxxxxxxxxxxxx
- 3: Gauze Pad: 2" x 2" (5 cm x 5 cm) origin:xxxxxxxxxxxxxxxxxxxx
- 1: Suture: 3-0 Braided Silk with Straight Needle origin:xxxxxxxxxxxxxxxxxxxx

¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU

California Prop. 65

WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex.

Store between 20 – 25°C (68 – 77°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



LOT LotNumber

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)LotNumber

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Set Product of xxx
Packaged in xxx



LBL059322 R03 (2024-05)



REF AK-04510-S

20 Ga. 12 cm catheter length .025 inch dia. spring-wire guide

Arterial Line

