



REF SAC-01618-CA

18<sup>Ga.</sup>16<sup>cm</sup>  
catheter  
length.025<sup>inch dia.</sup>  
spring-wire  
guide

# Arterial Catheterization

## Contents:

- |  |  |
|--|--|
| 1: Catheter: 18 Ga. x 16 cm Radiopaque PEBA  | 1: Suture: 3-0 Silk with Curved Needle                     |
| 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer | <sup>1</sup> A trademark of Becton, Dickinson and Company. |
| 1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW   | <sup>2</sup> A registered trademark of 3M Company.         |
| 1: Injection Needle: SafetyGlide™ <sup>1</sup> 25 Ga. x 1" (2.54 cm)   |  |
| 1: Injection Needle: SafetyGlide™ <sup>1</sup> 22 Ga. x 1-1/2" (3.81 cm) RW  |  |
| 1: Syringe: 3 mL Luer-Lock   |  |
| 1: Syringe: 5 mL Luer-Slip   |  |
| 1: SharpsAway® Disposal Cup  |  |
| 1: SharpsAway® II Locking Disposal Cup   |  |
| 1: Drape: 41" x 55" with 4" fenestration   |  |
| 1: Towel: 13" x 18" (33 cm x 46 cm)  |  |
| 1: Needle Holder   |  |
| 1: Safety Scalpel: #11   |  |
| 1: Checklist Sheet   |  |
| 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)  |  |
| 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)  |  |
| 1: Dressing: Tegaderm® <sup>2</sup> 10 cm x 12 cm  |  |
| 1: Surgical Apparel: Impervious Gown   |  |
| 1: Surgical Apparel: Bouffant Cap  |  |
| 1: Surgical Apparel: Mask with Eye Shield  |  |

---

**Catheter Priming Volume: 0.24 mL Priming volume is approximate.**

---

**Not made with natural rubber latex.**

